VULPES GOODS
BODYCARE

RED LIGHT THERAPY BELT

**MANUAL • ANLEITUG
HANDLEIDING • MANUEL**

WWW.VULPESGOODS.COM

ENGLISH	3
The Device	4
User Guide	6
Characteristics and Application	7
DEUTSCH	8
Das Gerät	9
Benutzerhandbuc	11
Eigenschaften und Anwendung	12
NEDERLANDS	13
Het apparaat	14
Gebruiksaanwijzing	16
Kenmerken en Toepassingen	17
FRANCAIS	18
Le dispositif	19
Manuel de l'utilisateur	20
Caractéristiques et applications	21

ENGLISH

THE DEVICE

Red light belt

Timer

Brightness + / -

On/Of

Band Switching

Frequency setting

Red Light Belt

Place the belt on your back, shoulders, or waist and adjust it to your preference. Ensure it is securely fastened with the Velcro closure for optimal contact and effect.

Band Switching

The band switching function allows you to switch between settings with each press of the button: first for ALL (850+660+850 nm), second for 660 nm red light, and third for 850 nm infrared light.

Timer

Press the timer button once to set the timer to 5 minutes. Press twice to set it to 10 minutes, three times for 15 minutes, four times for 20 minutes, five times for 25 minutes, and six times for 30 minutes.

Brightness + / -

The device has 5 brightness levels. Each press of the + / - button increases/decreases the brightness.

Frequency Setting

The frequency is adjustable from 10 Hz to 40 Hz, with each press of the button gradually increasing the frequency.

On/Off

Press the power button once to turn the device on. Press twice to turn the device off.

USER GUIDE

- 1.** Connect the Red Light belt and plug it into the power outlet.
- 2.** Press the power button once to turn the device on. Press again to turn the device off.
- 3.** Use the + / - button to select the desired brightness.
- 4.** Press the timer button to set the desired time.
- 5.** Press the frequency button to set the desired frequency (10Hz to 40Hz).
- 6.** Turn off the Red Light belt before unplugging it.

CHARACTERISTICS AND APPLICATION

- The light is not waterproof and should not be used in water or damp locations.
- Do not use the product with essential oils.
- To avoid the risk of electric shock, turn off the power switch before performing any installation.
- It is strictly forbidden to disassemble the device.
- Avoid using it on open wounds or irritated skin to prevent further damage or infection. Consult a medical specialist before use.
- Keep the device away from children to avoid accidents. It is not designed for use by children.
- While the device may be effective in relieving pain and promoting recovery, it is not a substitute for medical treatment. Consult a medical specialist for severe or persistent pain.

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